



## FlowSyn Polar Bear™

Advanced Cooling for Flow Chemistry

 **Advion Interchim**  
scientific

**UNIQSIS**  
*Accessible flow chemistry*

# FlowSyn Polar Bear™

Advanced cooling for flow chemistry



Clean, efficient cooling at the touch of a button

No heat transfer fluids, solid CO<sub>2</sub>, or solvents

Just plug in, switch on and set the temperature

FlowSyn Polar Bear™ is a state of the art chiller unit for low temperature flow-through chemistry applications.

**Powerful:** The reactor can be rapidly cooled down to any temperature between ambient and  $-88^{\circ}\text{C}$ .

**Rapid Cool-down:** Advanced cooling technology delivers much more rapid cooling than would be possible using a conventional chiller.


**Flexible:** The coil reactor module uses standard Uniqsis reactors and can accommodate single coil reactors up to 60 ml in volume, multiple smaller coils (available in a variety of sizes), and glass static mixer chips to deliver maximum flexibility for the flow chemist.

**Reagent Pre-cooling:** Reagent solutions can be efficiently pre-cooled prior to mixing by using either a static mixer chip or the dedicated heat exchange coils supplied with the interchangeable coil reactor top block.

**No 'Icing-up':** A vacuum-jacketed glass cover fitted with a nitrogen purge prevents ice formation and ensures that the reactor remains clearly visible at all times.

**Control via FlowSyn Interface:** FlowSyn Polar Bear™ can be controlled directly via the FlowSyn user interface. For example, in conjunction with FlowSyn Multi-X™ series of experiments, each at a different sub ambient temperature, can be performed in sequence.

UQ-1052 FlowSyn Polar Bear™ Specification	
Temperature range	ambient to $-88^{\circ}\text{C}$
Stability	$\pm 0.1^{\circ}\text{C}$ @ $-70^{\circ}\text{C}$ over 2 h
Power supply	220V 1100VA or 110V 1800VA
Dimensions	650 mm (w) x 650 mm (d) x 550 mm (h)
Weight	65 kg unpacked



Uniqsis Ltd, 29 Station Road, Shepreth, Cambridgeshire, SG8 6GB, United Kingdom  
+44 (0)845 864 7747 info@uniqsis.com www.uniqsis.com